

Examion

New

EXAMION DR810 PORTABLE

Portable DR Solution

*with mobile radiography
for the veterinary horse
practitioner*



V191400

EXAMION combi workstation accessories
including carrying case

SPECIFICATIONS

DETECTOR TYPE	FLAT-PANEL DETECTOR (DIRECT CONVERSION) TFT - amorphous selenium
DETECTING AREA	198 x 264 mm /8 x 10"
RESOLUTION	1,536 x 2,048 (3.1 M pixels)
PIXEL SIZE	129 x 129 µm
FILL FACTOR	~ 82 %
DYNAMIC RANGE	16,384 :1
A/D	14 bits
IMAGE ACQUISITION TIME	< 5 sec for full image
X-RAY ENERGY RANGE	50 ~ 150 kVp
X-RAY PULSE DURATION	1 ~ 20,000 ms
DIMENSIONS & WEIGHTS	PORTABLE (W X L X H): 331 x 377 x 22 mm, 3.0 kg approx. INDOOR (W X L X H): 369 x 381 x 32, 2.6 kg approx.
POWER CONSUMPTION	DC 12V 1.2A