

Examion


DR1417

New

EXAMION DR1417 RETROFIT

DR Retrofit Solution

*for radiological applications in the
veterinary practice*



Flat panel detector
344 x 429 mm

V191410

- ✓ **Direct Conversion Technology**
- ✓ Image viewing within seconds
- ✓ Very high image quality
- ✓ Incl. **EXAMION VETERINARY AQS SOFTWARE**

Examion

DR1417

New

SPECIFICATIONS

DETECTOR	TECHNOLOGY: TFT – amorphous selenium
	RESOLUTION: 2560 x 2048 pixels
	ACTIVE DETECTOR AREA: 344 x 429 mm / approx. 14 x 17"
	PIXEL SIZE: 168 µm
	A/D CONVERT: 14 bits
	DYNAMIC RANGE: 16384 : 1
	DATA ACQUISITION: < 3.6 sec. after x-ray image
	VOLTAGE: 100 V, 120 V, 230/240V (50/60 Hz)
	MAINS CONNECTION: 12 V, 1.5 A max.
	DETECTOR DIMENSIONS: 482 x 403 x 35 mm
	DETECTOR WEIGHT: 5.5 kg
WORKSTATION	EXAMION CS2 PREVIEW STATION PC with 19" touchscreen monitor



Includes the preview station:
PC with 19" touchscreen monitor