

Examion

New

EXAMION FLEXIBLE FPS

Wireless detector

Portable DR solution for use in the equine veterinary practice.


DR retrofit solution for use in the small animal veterinary practice.

GLOBAL INNOVATION



V191422

- ✓ Cassette format 35 x 43 cm
- ✓ Image viewing within seconds
- ✓ Optimal image quality
- ✓ Intuitive operation
- ✓ Incl. **EXAMION VETERINARY AQS SOFTWARE**



PORTABLE DR SOLUTION

DR RETROFIT SOLUTION

FLEXIBLE FPS

Compact, portable, wireless DR solution

- ✓ Wireless detector in cassette format 35 x 43 cm
- ✓ Rugged carrying case
- ✓ Intended for equine practices
- ✓ Combination workstation: Mini PC with 19" touchscreen monitor



FLEXIBLE FPS RETROFIT

Compact, portable, wireless DR solution

- ✓ Wireless detector in cassette format 35 x 43 cm
- ✓ Quick and easy DR retrofitting for stationary x-ray equipment
- ✓ Preview station: PC with 19" touchscreen monitor



SPECIFICATIONS

DETECTOR	<p>EXAMION FLEXIBLE FPS ISO 4090, Wireless Standard 802.11 n</p> <p>IMAGE SIZE: 35 x 43 cm</p> <p>MAXIMUM EXTERNAL DIMENSIONS: Corresponds to cassette size 35 x 43 cm</p> <p>WEIGHT: 3.9 kg</p> <p>PIXEL PITCH: 139 µm</p> <p>POWER SUPPLY: Integrated rechargeable battery using optional cable connection</p> <p>BATTERY LIFE: 70 to 90 images on one charge</p> <p>GRIDS: Can be used with standard grids</p>
BATTERY CHARGER	<p>POWER SUPPLY: 90-264 V, 50-60 Hz, max. 75W For simultaneous charging of up to three batteries incl. replacement batteries</p>
OPTIONAL POWER SUPPLY / COMMUNICATION CABLE	<p>For connecting the detector to the image viewer for image transmission per cable to the battery charger for trickle charging if the detector is used in a detector unit.</p>