

Examion

New

EXAMION MAXIVET 300HF

X-ray machine

Stationary CR/film x-ray machine
for radiological applications in the
small animal veterinary practice.



V191425

- ✓ 125 kV and 300 mA
- ✓ High frequency
- ✓ Floating table top, operated electrically with footswitch
- ✓ Capacitor technology: 3.5 kW input → 30 kW output!
- ✓ Socket-capable
- ✓ X-ray emitter electrically height-adjustable
- ✓ 605 organ programmes

SPECIFICATIONS

HIGH FREQUENCY X-RAY DEVICE	OUTPUT POWER: 30 kW
	VOLTAGE RANGE: 40 - 125 kV
	CURRENT RANGE: 1.0 - 100 mAs
	SWITCHING TIME: 0.001 - 1.25 s
	MAINS CONNECTION: 230 V - 50 Hz 1-phase (socket-capable)
	CAPACITOR TECHNOLOGY: 3.5 kW input → 30 kW output
	FOOT RELEASE SWITCH: 2-stage with safety function
	TABLE HEIGHT: 850 mm
	TABLE TOP: 144 x 66 cm, offset ± 38 cm
	TUBE: Electrically height-adjustable
TUBE	Toshiba x-ray tube with rotary anode and dual focus 1.0/2.0 mm
COLLIMATOR	Manual rotation, halogen lamp, timer
ACCESSORIES	Incl. scattered radiation grid 40 L/cm, 35 x 43 cm

